

# Nexus Medical Catalog

Catalog contains all Nexus products except neonatal/pediatric (see Nexus Neonatal/Pediatric Catalog)

Contact Nexus for pricing and availability; orders@nexusmedical.com

Bulk non-sterile available for most products, contact Nexus for details; orders@nexusmedical.com

Natural rubber latex is not used in the manufacturing of these products and tubing sets do not contain DEHP.



Sterile Part #	Description	Bonded	Pressure Rated	Rotating Male Luer	Clamp	APV	Qty/Box
NX4650	Nexus TKO <sup>®</sup> -6P Luer-Activated Device	N/A	325 psi	N/A	N/A	0.13 mL	100



N4685MP-B	6" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	325 psi	Yes	Pinch	0.4 mL	100
N4685MP	6" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	No	325 psi	Yes	Pinch	0.4 mL	100
N4685MN-B	6" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	325 psi	Yes	Slide	0.4 mL	100
N4685MN	6" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	No	325 psi	Yes	Slide	0.4 mL	100



N4685P-B	6" Standardbore Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	325 psi	Yes	Pinch	0.72 mL	100
N4685P	6" Standardbore Extension Set with Nexus TKO <sup>®</sup> -6P	No	325 psi	Yes	Pinch	0.72 mL	100
N4685N-B	6" Standardbore Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	325 psi	Yes	Slide	0.72 mL	100
N4685N	6" Standardbore Extension Set with Nexus TKO <sup>®</sup> -6P	No	325 psi	Yes	Slide	0.72 mL	100



N4695MN-B	6" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	325 psi	Yes	Slide	0.4 mL	100
N4695MN	6" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	No	325 psi	Yes	Slide	0.4 mL	100



N4683P-B	15" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	No	No	Pinch	0.45 mL	100
N4683P	15" Microbore Extension Set with Nexus TKO <sup>®</sup> -6P	No	No	No	Pinch	0.45 mL	100



N4690MN-B	6" Microbore Bifurcated Extension Set with (2) Nexus TKO <sup>®</sup> -6P	Yes	No	Yes	Slide	0.4 mL	100
N4690MN	6" Microbore Bifurcated Extension Set with (2) Nexus TKO <sup>®</sup> -6P	Yes	No	Yes	Slide	0.4 mL	100



N4691MN-B	6" Microbore Bifurcated Extension Set with (2) Nexus TKO <sup>®</sup> -6P and Power	Yes	325 psi	Yes	Slide	0.4 mL	100
N4691MN	6" Microbore Bifurcated Extension Set with (2) Nexus TKO <sup>®</sup> -6P and Power	No	325 psi	Yes	Slide	0.4 mL	100
N4691N-B	6" Standardbore Bifurcated Extension Set with (2) Nexus TKO <sup>®</sup> -6P and	Yes	325 psi	Yes	Slide	1.06 mL	100
N4691N	6" Standardbore Bifurcated Extension Set with (2) Nexus TKO <sup>®</sup> -6P and	No	325 psi	Yes	Slide	1.06 mL	100



NX4650HV	Nexus TKO <sup>®</sup> -6PHV Hi-Visibility Luer-Activated Device	N/A	325 psi	N/A	N/A	0.13 mL	100
----------	--	-----	---------	-----	-----	---------	-----



N4685MPHV	7" Microbore Extension Set with Nexus TKO <sup>®</sup> -6PHV Hi-Visibility	No	No	Yes	Pinch	0.32 mL	100
N4683PHV-B	15" Microbore Extension Set with Nexus TKO <sup>®</sup> -6PHV Hi-Visibility	Yes	No	No	Pinch	0.45 mL	100
N4683PHV	15" Microbore Extension Set with Nexus TKO <sup>®</sup> -6PHV Hi-Visibility	No	No	No	Pinch	0.45 mL	100



# Nexus Medical Catalog

Catalog contains all Nexus products except neonatal/pediatric (see Nexus Neonatal/Pediatric Catalog)

Contact Nexus for pricing and availability; orders@nexusmedical.com

Bulk non-sterile available for most products, contact Nexus for details; orders@nexusmedical.com

Natural rubber latex is not used in the manufacturing of these products and tubing sets do not contain DEHP.



Sterile Part #	Description	Pressure Rated	Rotating Male Luer	Clamp	APV	Qty/Box
N3560	8" 2-up Manifold with Nexus Two-Way Valves	No	Yes	Pinch	0.1 mL (ea valve) 0.9 mL (set)	100
N3560-1W	8" 2-up Manifold with Nexus One-Way Check Valves	No	Yes	Pinch	0.07 mL (ea valve) 0.9 mL (set)	100



N3561	10" 3-up Manifold with Nexus Two Way Valves	No	Yes	Pinch	0.1 mL (ea valve) 1.0 mL (set)	100
N3561-1W	10" 3-up Manifold with Nexus One-Way Check Valves	No	Yes	Pinch	0.07 mL (ea valve) 1.0 mL (set)	100



N3574	10" 4-up Manifold with Nexus Two-Way Valves	No	Yes	Pinch	0.1 mL (ea valve) 0.9 mL (set)	50
N3574-1W	10" 4-up Manifold with Nexus One-Way Check Valves	No	Yes	Pinch	0.07 mL (ea valve) 0.9 mL (set)	50



N3576	10" 6-up Manifold with Nexus Two-Way Valves	No	Yes	Pinch	0.1 mL (ea valve) 0.9 mL (set)	50
N3576-1W	10" 6-up Manifold with Nexus One-Way Check Valves	No	Yes	Pinch	0.07 mL (ea valve) 0.9 mL (set)	50




# Nexus Medical Catalog

Catalog contains all Nexus products except neonatal/pediatric (see Nexus Neonatal/Pediatric Catalog)

Contact Nexus for pricing and availability; orders@nexusmedical.com

Bulk non-sterile available for most products, contact Nexus for details; orders@nexusmedical.com

Natural rubber latex is not used in the manufacturing of these products and tubing sets do not contain DEHP.



Sterile Part #	Description	Bonded	Pressure Rated	Rotating Male Luer	Clamp	APV	Qty/Box
N6500	Nexus NIS®-6P Luer-Activated Device	N/A	325 psi	N/A	N/A	0.08 mL	100



N6966MP-B	6" Microbore Extension Set w/ Nexus NIS®-6P	Yes	325 psi	Yes	Pinch	0.3 mL	100
N6966MP	6" Microbore Extension Set w/ Nexus NIS®-6P	No	325 psi	Yes	Pinch	0.3 mL	100
N6966MN-B	6" Microbore Extension Set w/ Nexus NIS®-6P	Yes	325 psi	Yes	Slide	0.3 mL	100
N6966MN	6" Microbore Extension Set w/ Nexus NIS®-6P	No	325 psi	Yes	Slide	0.3 mL	100



N6966P-B	6" Standardbore Extension Set w/ Nexus NIS®-6P	Yes	325 psi	Yes	Pinch	0.65 mL	100
N6966P	6" Standardbore Extension Set w/ Nexus NIS®-6P	No	325 psi	Yes	Pinch	0.65 mL	100
N6966N-B	6" Standardbore Extension Set w/ Nexus NIS®-6P	Yes	325 psi	Yes	Slide	0.65 mL	100
N6966N	6" Standardbore Extension Set w/ Nexus NIS®-6P	No	325 psi	Yes	Slide	0.65 mL	100



N6980N-B	8" Minibore Extension Set w/ Nexus NIS®-6P	Yes	No	No	Slide (removable)	0.38 mL	100
N6980N	8" Minibore Extension Set w/ Nexus NIS®-6P	No	No	No	Slide (removable)	0.38 mL	100



N6983P-B	15" Microbore Extension Set w/ Nexus NIS®-6P	Yes	No	No	Pinch	0.4 mL	100
N6983P	15" Microbore Extension Set w/ Nexus NIS®-6P	No	No	No	Pinch	0.4 mL	100



N6976MN-B	7" Microbore Bifurcated Extension Set w/ (2) Nexus NIS®-6P	Yes	No	Yes	Slide	0.34 mL	100
N6976MN	7" Microbore Bifurcated Extension Set w/ (2) Nexus NIS®-6P	No	No	Yes	Slide	each leg	100
N6976N-B	7" Standardbore Bifurcated Extension Set w/ (2) Nexus NIS®-6P	Yes	No	Yes	Slide	1.0 mL	100
N6976N	7" Standardbore Bifurcated Extension Set w/ (2) Nexus NIS®-6P	No	No	Yes	Slide	each leg	100



# Nexus Medical Catalog

Catalog contains all Nexus products except neonatal/pediatric (see Nexus Neonatal/Pediatric Catalog)

Contact Nexus for pricing and availability; orders@nexusmedical.com

Bulk non-sterile available for most products, contact Nexus for details; orders@nexusmedical.com

Natural rubber latex is not used in the manufacturing of these products and tubing sets do not contain DEHP.



Sterile Part #	Description	Bonded	Pressure Rated	Rotating Male Luer	Clamp	APV	Qty/Box
N8000	Nexus NIS <sup>®</sup> -8 Split Septum Device	N/A	325 psi	N/A	N/A	0.1 mL	100



N8046	4" Minibore Extension Set w/ Split Septum T-Site	N/A	No	Yes	Pinch	0.25 mL	100
N8066	6" Minibore Extension Set w/ Split Septum T-Site	N/A	No	Yes	Pinch	0.3 mL	100



P4842-B	4" Split Septum T-Site Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	No	Yes	Pinch	0.36 mL	100
P4842	4" Split Septum T-Site Extension Set with Nexus TKO <sup>®</sup> -6P	No	No	Yes	Pinch	0.36 mL	100
P4862-B	6" Split Septum T-Site Extension Set with Nexus TKO <sup>®</sup> -6P	Yes	No	Yes	Pinch	0.48 mL	100
P4862	6" Split Septum T-Site Extension Set with Nexus TKO <sup>®</sup> -6P	No	No	Yes	Pinch	0.48 mL	100



N0656	Lever Lock Cannula	N/A					200
N0657	Straight Blunt Cannula (blister pack)	N/A					100



# Nexus Medical Catalog

Catalog contains all Nexus products except neonatal/pediatric (see Nexus Neonatal/Pediatric Catalog)

Contact Nexus for pricing and availability; orders@nexusmedical.com

Bulk non-sterile available for most products, contact Nexus for details; orders@nexusmedical.com

Natural rubber latex is not used in the manufacturing of these products and tubing sets do not contain DEHP.



Sterile Part #	Description	Bonded	Pressure Rated	Rotating Male Luer	Clamp	APV	Qty/Box
N0002	7" Microbore Extension Set	N/A	No	No	Slide (removable)	0.2 mL	100
N0003	12" Microbore Extension Set	N/A	No	No	Slide (removable)	0.3 mL	100
N0006	36" Microbore Extension Set	N/A	No	No	Slide (removable)	0.8 mL	100
N0007	60" Microbore Extension Set	N/A	No	No	Slide (removable)	1.2 mL	100
N1602	60" Minibore Extension Set	N/A	No	No	Slide (removable)	2.0 mL	100
N1151	15" Microbore Extension Set	N/A	No	No	Pinch	0.4 mL	100
N0008	60" Ultra-Microbore Extension Set	N/A	No	No	No Clamp	0.4 mL	100
N1602A	60" Minibore Extension Set	N/A	No	No	Slide	2.0 mL	100



N5108	10" Standardbore Ext Set w/ SmartSite® Y-Connector	N/A	No	Yes	Slide	1.2 mL	50
-------	--	-----	----	-----	-------	--------	----



N1516	35" Standardbore Set with (2) Nexus Y-Connectors	N/A	No	Yes	Slide	5.5 mL	100
-------	--	-----	----	-----	-------	--------	-----



N1081	7" Microbore Extension Set with .22 Micron Filter	N/A	No	No	Slide	0.44 mL	100
N1082	60" Ultra-Microbore Extension Set with .22 Micron Filter	N/A	No	No	Slide	0.65 mL	100
N1083	60" Minibore Extension Set with 1.2 Micron Filter	N/A	No	No	Slide	2.3 mL	100



filters can be round or rectangle on actual sets

N0080	12" Standardbore Extension Set with .22 Micron Filter	N/A	No	Yes	Slide	3.8 mL	50
N0086	12" Standardbore Extension Set with 1.2 Micron Filter	N/A	No	Yes	Slide	3.8 mL	50



BD9001	Single Entry Trifurcated Blood Sampling Set	N/A	No	Yes	Pinch	0.16 mL (each leg)	100
--------	---	-----	----	-----	-------	--------------------	-----



# Nexus Medical Catalog

Catalog contains all Nexus products except neonatal/pediatric (see Nexus Neonatal/Pediatric Catalog)

Contact Nexus for pricing and availability; orders@nexusmedical.com

Bulk non-sterile available for most products, contact Nexus for details; orders@nexusmedical.com

Natural rubber latex is not used in the manufacturing of these products and tubing sets do not contain DEHP.



Sterile Part #	Description	Bonded	Pressure Rated	Rotating Male Luer	Clamp	APV	Qty/Box
N9902R	Red End Cap Male/Female			N/A			see note below
N9902B	Blue End Cap Male/Female						



**sold as CASE ONLY (36 boxes of 200 per case)**

N0650	20mm Multi-Dose Vial Access Device	Yes	No		N/A		50
-------	------------------------------------	-----	----	--	-----	--	----



N0670	Non-Vented Bag Spike with Nexus NIS <sup>®</sup> -6P	Yes	No		N/A		100
-------	--	-----	----	--	-----	--	-----



N0606	Mini Vented Transfer Pin				N/A		50
-------	--------------------------	--	--	--	-----	--	----



BD6010	Blood Collection Assembly with Male Luer Lock				N/A		200
BD6018	Blood Transfer Assembly with Female Luer Lock						200



# Nexus Medical Catalog

Catalog contains all Nexus products except neonatal/pediatric (see Nexus Neonatal/Pediatric Catalog)

Contact Nexus for pricing and availability; orders@nexusmedical.com

Bulk non-sterile available for most products, contact Nexus for details; orders@nexusmedical.com


Natural rubber latex is not used in the manufacturing of these products and tubing sets do not contain DEHP.




Sterile Part #	Description	Bonded	Pressure Rated	Rotating Male Luer	Clamp	APV	Qty/Box
N1512	96" I.V. Administration Set with (2) SmartSite® Needle Free Y-Sites, 15 Drop m/L Universal Drip Chamber, Check Valve	N/A	No	Yes	Roller/Slide	18 mL	50




N1523	92" I.V. Administration Set with Flow Controller, 15 drops/mL Universal Drip Chamber, 1 Injection Y-Site,	N/A	No	Yes	Pinch	18 mL	50
-------	---	-----	----	-----	-------	-------	----



N1540A	40" I.V. Secondary Set w/ Hanger, 15 drops/mL Universal Drip Chamber	N/A	No	Yes	Roller	8 mL	50
--------	--	-----	----	-----	--------	------	----



N1545	109" I.V. Administration Set with (3) SmartSite® Needle-Free Y-Sites, 15 drops/mL Universal Drip Chamber, Check Valve	N/A	No	Yes	Roller/Slide	20 mL	50
-------	---	-----	----	-----	--------------	-------	----



N2523	20" I.V. Extension Set with Flow Controller, Injection Y-Site	N/A	No	Yes	N/A	3 mL	50
-------	---	-----	----	-----	-----	------	----

